




Imaging



Molecular Imager® Gel Documentation Systems



The Worldwide Leader in Gel Imaging





Green Fluorescent Protein



Silver stain



Immun-Star™ HRP kit

The Unmistakable Leader in High-Resolution Systems

Bio-Rad is the worldwide leader when it comes to providing high-quality imaging products. Whether your research includes everyday imaging of PCR products, more demanding low-light applications such as chemiluminescent western blots, or both — the Gel Doc™ XR and ChemiDoc™ XRS imaging systems meet your needs.

Two Complete Gel Documentation Systems

The Gel Doc XR and ChemiDoc XRS systems ship complete with darkroom, camera, UV and white light illumination, filter slider with amber filter, UV protection shield, Quantity One® 1-D analysis software, and unlimited copies of Quantity One® Basic software for data sharing.



Extra
Resolution

Gel Doc XR

The Gel Doc XR system is a high-resolution gel documentation system for all of your everyday imaging needs.

The 1.4 megapixel imaging system is fast and easy to use. Follow the onscreen prompts and get your results with only three clicks of a mouse. The motorized zoom lens allows a hands-free approach to finding the perfect focus, zoom, and iris settings. The intuitive and easy-to-use Quantity One interface allows you to expose and print directly from the acquisition window. The Gel Doc XR provides a detection range of 3 orders of magnitude and utilizes a FireWire interface for fast data transfer.

- High-resolution imaging with 1.4 million pixels
- Motorized zoom lens for hands-free gel documentation
- Real-time imaging for quick positioning and focusing of samples
- Onscreen integration — no need to save an image; simply “freeze” the image at the desired intensity and print
- FireWire interface for rapid data transfer
- White light transilluminator allows documentation of visible dyes
- Upgradable to the ChemiDoc XRS system

The Gel Doc XR provides extra resolution for the acquisition, optimization, quantitation, and documentation of all of your gel data.



Extra Resolution
and Sensitivity

ChemiDoc XRS

The ChemiDoc XRS system is an advanced chemiluminescent detection system offering extra resolution and sensitivity.

The system features flat fielding technology for superior uniformity and quantitation, and the flexibility to image chemiluminescent, fluorescent, and colorimetric samples. The system includes a supersensitive CCD, supercooled for excellent detection of faint samples. The ChemiDoc XRS can be used for imaging in a wide variety of applications that require high resolution (such as 2-D gel electrophoresis) and excellent sensitivity (such as chemiluminescent western blots).


- Supercooled high-resolution CCD camera with 1.44 megapixel resolution
- User control of camera bin settings
- Exceptional dynamic range and sensitivity
- Real-time imaging
- Dynamic flat fielding
- Excellent collection efficiency — f/1.2 motorized zoom lens or f/0.95 wide-angle lens

The ChemiDoc XRS is optimized for high-resolution chemiluminescent and fluorescent imaging.

SYPRO Ruby stain

SYBR Green I stain

Coomassie Blue R-250 stain



Quantity One 1-D Analysis Software

Quantity One software is a powerful package for imaging and analyzing 1-D electrophoretic gels and blots, dot blots, slot blots, and colonies. Quantity One software analyzes a variety of data collected from densitometers, storage phosphor imagers, fluorescence imagers, and gel documentation systems. Quantity One software automates analysis for fast, high-quality results.

User-Friendly Interface

- Built-in Quick Guides
- Large selection of predetermined settings for the most commonly used detection reagents
- Annotation tools

Analysis Tools

- Automatic band and lane detection
- Molecular weight determination
- Colony counting and volume analysis
- Automation Manager for archiving and retrieving your most commonly used analysis templates
- 3-D viewer for critical analysis of closely spaced bands

Data Output

- Data sharing with unlimited copies of Quantity One Basic
- TIFF file import, export, and publishing
- JPEG file export
- Mac and PC compatibility
- Networkability

Now With Data Security

- Tools for compliance with US FDA 21 CFR Part 11 regulations



Molecular Imager **Gel Documentation Systems**

Specifications

	Gel Doc XR System	ChemiDoc XRS System
CCD resolution (H x V)	1,360 x 1,024 μm	1,392 x 1,040 μm
Pixel density	4,096 gray levels	65,536 gray levels
Pixel size (H x V)	4.6 x 4.6 μm	6.7 x 6.7 μm
Resolution	1.4 megapixel	1.44 megapixel
Dynamic range	>3 orders of magnitude	4 orders of magnitude
Camera cooling system	—	Peltier
Camera cooling temperature	—	-50°C
Motorized zoom lens	C-mount, f/1.2, 8.5–51 mm	C-mount, f/1.2, 12–75 mm
Illumination modes	Trans-UV, white, epi-white	Trans-UV, white, epi-white
Excitation source	254, 302, 365 nm, and white light	254, 302, 365 nm, and white light
Filter positions	Fluorescence: 2	Fluorescence: 2; chemiluminescence: 1
Emission filters	1 included (amber), 4 optional	1 included (amber), 4 optional
Transillumination area	25 x 26 cm	25 x 26 cm
Dynamic flat fielding	—	Yes, CV* \leq 5%
Operating system compatibility	Windows and Mac OS	Windows and Mac OS
Dimensions (W x D x H)	60 x 36 x 96 cm	60 x 36 x 96 cm
Weight	32 kg (70.4 lb)	32 kg (70.4 lb)

* Coefficient of variance.

Ordering Information

Catalog #	Description
Gel Documentation Systems	
170-8170	Molecular Imager Gel Doc XR System, PC
170-8171	Molecular Imager Gel Doc XR System, Mac
170-8070	Molecular Imager ChemiDoc XRS System, PC
170-8071	Molecular Imager ChemiDoc XRS System, Mac
170-9600	Quantity One 1-D Analysis Software
170-9615	Quantity One CFR Module

Accessories

170-7950	White Light Transilluminator, universal hood
170-8001	White Light Conversion Screen
170-8072*	f/0.95 25 mm Wide-Angle Lens
170-8073*	f/0.95 17 mm Wide-Angle Lens
170-8074	Filter, 520DF30 62 mm for SYBR Green/GFP/SYBR Gold/fluorescein
170-8075	Filter, 560DF50 62 mm for Cy3/rhodamine
170-8076	Filter, 630BP30 62 mm for SYPRO Ruby/Texas Red
170-8077	Filter, 480BP70 62 mm for Hoechst/coumarin
170-8098	254 nm UV Lamps, 6
170-6887	365 nm UV Lamps, 6
170-8097	Standard 302 nm UV Lamps, 6
170-8089	Mitsubishi P93DW Printer
170-7581	Mitsubishi Thermal Printer Paper, 4 rolls
170-8067	Sony UPD895 Printer Paper, 5 rolls

Computers

Inquire

* For ChemiDoc XRS system only.

Coomassie is a trademark of BASF Aktiengesellschaft. Cy is a trademark of GE Healthcare. FireWire, Mac, and Mac OS are trademarks of Apple Computer, Inc. Hoechst is a trademark of Hoechst Aktiengesellschaft. Mitsubishi is a trademark of Mitsubishi Electric Corporation. Sony is a trademark of Sony Corporation. SYBR, SYPRO, and Texas Red are trademarks of Molecular Probes, Inc. Windows is a trademark of Microsoft Corporation. Practice of the polymerase chain reaction (PCR) may require a license.



BIO-RAD

**Bio-Rad
Laboratories, Inc.**

Life Science
Group

Web site www.bio-rad.com **USA** 800 4BIORAD **Australia** 02 9914 2800 **Austria** 01 877 89 01 **Belgium** 09 385 55 11 **Brazil** 55 21 3237 9400
Canada 905 712 2771 **China** 86 21 6426 0808 **Czech Republic** 420 241 430 532 **Denmark** 44 52 10 00 **Finland** 09 804 22 00 **France** 01 47 95 69 65
Germany 089 318 84 0 **Greece** 30 210 777 4396 **Hong Kong** 852 2789 3300 **Hungary** 36 1 455 8800 **India** 91 124 4029300/5013478 **Israel** 03 963 6050
Italy 39 02 216091 **Japan** 03 5811 6270 **Korea** 82 2 3473 4460 **Mexico** 55 5200 05 20 **The Netherlands** 0318 540666 **New Zealand** 64 9415 2280
Norway 23 38 41 30 **Poland** 48 22 331 99 99 **Portugal** 351 21 472 7700 **Russia** 7 095 721 14 04 **Singapore** 65 6415 3188 **South Africa** 27 0861 246 723
Spain 34 91 590 5200 **Sweden** 08 555 12700 **Switzerland** 061 717 95 55 **Taiwan** 886 2 2578 7189/2578 7241 **United Kingdom** 020 8328 2000