






Experion[™]

Automated Electrophoresis System

Focus on the results, not the method.

BIO-RAD

Experience Meets Innovation



The Experion automated electrophoresis system is a powerful and affordable separation and analysis system that applies microfluidic technology to reinvent the way that you perform protein and RNA electrophoresis. The Experion system combines Bio-Rad's expertise in electrophoresis with the innovation of Caliper Life Sciences' LabChip technology to deliver new levels of performance in automation. The Experion system advances automated electrophoresis to expand your ability to produce data quickly, without compromising the quality of results.

Rapid, Automated Results

The Experion system automatically performs the multiple steps of gel-based electrophoresis. You can walk away and do more with your time while the Experion system produces highly reproducible separation and quantitation of your protein and RNA samples.

Integrated System Design

The sleek components of the Experion system incorporate efficient and creative designs to deliver high-quality results. Optimized microfluidic chip design, electrophoresis-grade reagents, exclusive protein and RNA standards, easy-to-use automated

priming and electrophoresis stations, and powerful software analysis tools all combine to form an integrated system that streamlines separation and analysis.

Major Benefits

- Dramatically reduced time-to-results, hands-on time, reagent usage, and sample consumption
- An affordable alternative to traditional electrophoresis
- The highest protein resolution and sensitivity and best quantitation results available in an automated system
- Accurate RNA quantitation

Automated Electrophoresis

Separation

Staining

Destaining

Imaging

Band detection

Quantitation

Basic data analysis

Results in only
30 minutes



System Features

Superior Performance

- Fast, 30 minute batch runs of protein and RNA samples
- Accurate single-step protein sizing from 10 to 260 kD
- 2-in-1 process for RNA: integrity checks and quantitation
- Exclusive protein and RNA standards produce accurate and reproducible sizing and quantitation
- Flexible software tools make data analysis easier and more efficient

Convenient Data Analysis Tools

- Sizing and quantitation calculations performed automatically
- Intuitive navigation of separation and data analysis screens
- Quick comparisons of protein or RNA components across the chip
- Regulatory features — tools for US FDA 21 CFR Part 11 compliance and installation qualification/operational qualification (IQ/OQ) functions

User-Friendly

- Automated and integrated system makes electrophoresis easier than ever before
- Automatic and error-free chip priming
- Minimal hands-on time required for unattended operation
- Minimal sample and reagent requirements
- Reduced exposure to hazardous chemicals



Automate Your Protein and RNA Analysis



Superior Analysis Kits and Chips Provide Improved Resolution and Quantitation

Experion Analysis Kits

Experion analysis kits combine state-of-the-art chip design with high-quality reagents to perform reproducible, quantitative, and accurate protein and RNA analysis in minutes. Streamlined chip preparation methods and minimal sample requirements result in rapid experiments with minimal hands-on time.

Each analysis kit includes:



- Experion Pro260, RNA StdSens, or RNA HighSens chips
- High-quality gel matrix for separation and resolution similar to mini gels
- High-sensitivity fluorescent dye for accurate detection
- Experion protein or RNA ladder for accurate sizing and quantitation
- Optimized sample buffer for accurate quantitation and reproducible results

Kits are available in flexible ordering configurations to match your research needs.

Experion Pro260 Analysis Kit

The Experion Pro260 analysis kit delivers fast, sensitive, and reproducible analyses of protein samples.

- Analysis of up to 10 samples in 30 minutes
- Resolution and quantitation of 10–260 kD proteins
- Improved resolution over other automated systems
- Sensitivity comparable to that of colloidal Coomassie Blue gel staining
- Protein sizing, quantitation, and analysis in a single step



Run time, sec



Comparison of separation of proteins from equal amounts of *E. coli* lysate on automated electrophoresis systems.

Upper panel, separation on a Pro260 chip, displayed using Experion software; lower panel, separation on a competitor's chip, displayed using a competitor's automated electrophoresis system. Run time, ~60 sec. Note the greater number of peaks and increased resolution of the Experion Pro260 result.

Experion RNA HighSens and RNA StdSens Analysis Kits

The Experion RNA analysis kits offer fast, accurate, and reproducible purity checks designed to provide confidence in RNA quality without delaying your experiments.

- Analysis of 1–12 samples in 10–30 minutes, depending on kit
- Quantitation at nanogram (RNA StdSens kit) amounts
- RNA ladder included in each kit
- Single-step RNA purity assessment



Comparison of RNA separation on automated systems.

Upper panel, Experion RNA ladder separated on an RNA StdSens chip, displayed using Experion software; lower panel, a competitor's RNA ladder separated and displayed using a competitor's automated electrophoresis system. Run time, ~60 sec. The Experion RNA ladder provides more uniform peak heights (fluorescence intensities), resulting in improved quantitation.



Experion Automated Electrophoresis Station

The Experion automated electrophoresis station performs all the steps of gel-based electrophoresis in one compact, durable unit. Its multifunctionality combines electrophoresis, staining, destaining, band detection, and imaging into a single 30 minute step.

Electrode manifold with 16 high-quality platinum pins for reproducible runs

Easy-to-access platform for inserting or removing chips



- Highly accurate laser provides precise fluorescence detection
- USB port allows easy installation and maximum connectivity
- Built-in power supply reduces cost and saves benchspace
- Large LED "on" light blinks to indicate a run is in progress



Software

Experion Priming Station

The Experion automated priming station consistently prepares protein and RNA chips for successful electrophoresis with minimal hands-on time. Preset time and pressure settings ensure optimal introduction, or priming, of the gel matrix into the microchannels of the chip.



- Large LCD display clearly shows the preset time and pressure settings
- Integrated timer conveniently counts down the time-sensitive priming step
- Accessible chip platform allows easy chip placement and sample loading
- Alignment arrows on chip and priming station ensure proper chip placement for successful priming
- Secure locking mechanism prevents early release while priming
- Built-in, pressure-activated release mechanism ensures precise priming


Experion Vortex Station

The Experion vortex station ensures complete mixing of RNA samples and analysis reagents for effective sample runs.



- Prongs on the vortex adaptor securely hold the chip in place
- Preset speed and time settings provide single-step, precise mixing of samples and reagents
- Mixing within the chip reduces reagent volume and pipetting steps
- Beveled edges on the vortex adaptor provide easy access while loading and unloading a chip

Tree view Electropherogram view Gel view




Experion software makes sophisticated technology surprisingly user-friendly. Experimental results are analyzed through easy-to-use menus and result tables.

Result table

Experion software is your entryway into automated electrophoresis. The simple yet comprehensive working screen and built-in analysis functions allow you to obtain the information you need without spending a lot of your valuable time.

Versatile and Efficient

- Perform the run and analyze the data from a single screen
- Electropherogram (peak) or gel views offer easy access to information in both formats
- Data are organized in a tree-view format for logical storage, sorting, and retrieval of run information




Experion Pro260 Chip

The Pro260 chip allows analysis of up to 10 protein samples (10–260 kD) in approximately 30 minutes.

Experion RNA HighSens Chip

The RNA HighSens chip allows analysis of up to 11 RNA samples in the 100–5,000 pg/µl range.



Experion software automatically calculates concentration, protein molecular weight, percent of total sample, and, in total RNA runs, the ribosomal RNA ratio. A typical total RNA analysis is shown above.


Experion software compares data from separate wells within a chip. The same protein peak from different wells is shown to the right.

W	Well ID	Sample Name	Mol. Wt. (kDa)	Concentration (µg/ml)	% Total Expression
1	1	E. coli (2 mg/ml)	43.66	241.61	7.50
2	2	E. coli (2 mg/ml)	43.34	218.67	6.98
3	3	E. coli (2 mg/ml)	43.60	202.50	6.12
4	4	E. coli (2 mg/ml)	43.85	217.88	6.72
5	5	E. coli (2 mg/ml)	43.55	215.56	6.39
6	6	E. coli (2 mg/ml)	43.49	211.43	6.24
7	7	E. coli (2 mg/ml)	43.85	182.79	5.68
8	8	E. coli (2 mg/ml)	44.02	192.79	5.47
9	9	E. coli (2 mg/ml)	43.87	203.48	5.80
Mean			43.85	209.42	6.34
Std. Deviation			0.19	16.24	0.60
% CV			0.44	7.76	9.46

Accessible Information

Experion software provides innovative tools that take the tedium out of data analysis.

- Automatic calculation of size, concentration, and percent of total sample — results for each peak appear in the result table
- Query-based comparisons of a single peak across all samples in a chip enable statistical analysis of the expression of a single protein or RNA of interest
- Data export function allows data sets to be exported to a spreadsheet for customized analysis



Experion RNA StdSens Chip

The RNA StdSens chip allows analysis of up to 12 RNA samples in the 5–500 ng/µl range.

Experion LabChip Technology

The Experion analysis chip houses LabChip technology developed by Caliper Life Sciences, Inc. The chip is a powerful, miniaturized device — much like a tiny laboratory — that combines the functionality of several larger benchtop analytical instruments. Up to 10 protein or 12 RNA samples can be analyzed in only 30 minutes.

Research Applications

The Experion automated electrophoresis system is the perfect complement to Bio-Rad's protein separation and gene expression analysis tools.

Protein Analysis

A wide variety of protein-related applications are supported by the Experion system, including quality control, protein purity and stability analysis, protocol optimization, and evaluation of recombinant protein expression.

Proteins of interest are often isolated and purified by fractionating a complex sample using chromatography systems, such as the BioLogic DuoFlow™ system. The Experion system is ideal for assessing any purification protocol.



RNA Analysis

Gene expression profiling experiments require high-purity RNA to ensure optimal results. The Experion system is perfect for evaluation of RNA purity because it requires very little of your valuable samples and your valuable time. The VersArray® microarray systems offer instruments for high-precision microarray experiments. As target transcripts are identified, the iQ™5 and MyiQ™ real-time PCR systems provide accurate transcript quantitation.



iQ5 Multicolor Real-Time PCR Detection System



MyiQ Single-Color Real-Time PCR Detection System

Ordering Information

Catalog #	Description
Experion Automated Electrophoresis Systems	
700-7000	Experion System , 100–120/220–240 V, for protein analysis, includes electrophoresis station, priming station, software, USB2 cable, instructions
700-7001	Experion System , 100–120 V, for RNA analysis, includes electrophoresis station, priming station, vortex station, software, USB2 cable, instructions
700-7002	Experion System , 220–240 V, for RNA analysis
Electrophoresis Station and Replacement Parts	
700-7010	Experion Electrophoresis Station , 100–120/220–240 V, includes USB2 cable, instructions
700-7020	Experion Electrode Manifold , replacement
700-7021	Experion Lid , replacement
700-7022	Experion USB2 Cable With Ferrite , replacement
Priming Station and Replacement Parts	
700-7030	Experion Priming Station , 100–120/220–240 V, includes 2 priming seals
700-7031	Experion Priming Seals , replacement, provides air seal on top of priming well, 2
Vortex Station and Replacement Parts (for RNA Analysis)	
700-7040	Experion Vortex Station , 115 V
700-7041	Experion Vortex Station , 230 V
700-7042	Experion Vortex Adaptor , holds analysis chip in vortex station, replacement
Experion Analysis Kits	
700-7101	Experion Pro260 Analysis Kit for 10 Chips , includes 10 Pro260 chips, 3 x 520 µl Pro260 gel, 45 µl Pro260 stain, 60 µl Pro260 ladder (10–260 kD), 400 µl Pro260 sample buffer, 3 spin filters
700-7102	Experion Pro260 Analysis Kit for 25 Chips , includes 25 Pro260 chips, 5 x 520 µl Pro260 gel, 2 x 45 µl Pro260 stain, 2 x 60 µl Pro260 ladder (10–260 kD), 2 x 400 µl Pro260 sample buffer, 5 spin filters
700-7103	Experion RNA StdSens Analysis Kit for 10 Chips , includes 10 RNA StdSens chips, 1,250 µl RNA gel, 20 µl RNA StdSens stain, 20 µl RNA ladder, 900 µl RNA StdSens loading buffer, 2 spin filters
700-7104	Experion RNA StdSens Analysis Kit for 25 Chips , includes 25 RNA StdSens chips, 2 x 1,250 µl RNA gel, 2 x 20 µl RNA StdSens stain, 2 x 20 µl RNA ladder, 2 x 900 µl RNA StdSens loading buffer, 4 spin filters
700-7105	Experion RNA HighSens Analysis Kit for 10 Chips , includes 10 RNA HighSens chips, 1,250 µl RNA gel, 20 µl RNA HighSens stain, 20 µl RNA ladder, 900 µl RNA HighSens loading buffer, 100 µl RNA sensitivity enhancer, 2 spin filters
700-7106	Experion RNA HighSens Analysis Kit for 25 Chips , includes 25 RNA HighSens chips, 2 x 1,250 µl RNA gel, 2 x 20 µl RNA HighSens stain, 20 µl RNA ladder, 2 x 900 µl RNA HighSens loading buffer, 2 x 100 µl RNA sensitivity enhancer, 4 spin filters
Experion Analysis Kit Accessories	
700-7151	Experion Pro260 Chips , 10
700-7152	Experion Pro260 Reagents and Supplies , for 10 chips, includes 3 x 520 µl Pro260 gel, 45 µl Pro260 stain, 60 µl Pro260 ladder (10–260 kD), 400 µl Pro260 sample buffer, 3 spin filters
700-7256	Experion Pro260 Ladder , 60 µl (10–260 kD)

Catalog #	Description
Experion Analysis Kit Accessories (Cont.)	
700-7153	Experion RNA StdSens Chips , 10
700-7154	Experion RNA StdSens Reagents and Supplies , for 10 chips, includes 1,250 µl RNA gel, 20 µl RNA StdSens stain, 20 µl RNA ladder, 900 µl RNA StdSens loading buffer, 2 spin filters
700-7155	Experion RNA HighSens Chips , 10
700-7156	Experion RNA HighSens Reagents and Supplies , for 10 chips, includes 1,250 µl RNA gel, 20 µl RNA HighSens stain, 20 µl RNA ladder, 900 µl RNA HighSens loading buffer, 100 µl RNA sensitivity enhancer, 2 spin filters
700-7255	Experion RNA Ladder , 20 µl
700-7251	Experion Cleaning Chips , 10
700-7252	Experion Electrode Cleaner , 250 ml
700-7253	Experion DEPC-Treated Water , 100 ml
700-7254	Experion Spin Filters , 10
Experion Software*	
700-7050	Experion Software , system operation and data analysis tools, PC
700-7051	Experion Validation Kit , 3 test chips, qualification procedures, dongle, PC
700-7052	Experion Security Edition Software , system operation, standard and 21 CFR Part 11 data analysis tools, 3 test chips, qualification procedures, dongle, PC

* Optional computer systems available. Contact your local Bio-Rad representative for more information, including specific computer requirements. Or visit us on the Web at www.bio-rad.com/experion/



LabChip and the LabChip logo are trademarks of Caliper Life Sciences, Inc. Bio-Rad Laboratories, Inc. is licensed by Caliper Life Sciences, Inc. to sell products using the LabChip technology for research use only.

These products are licensed under US Patent Nos. 5,863,753, 5,658,751, 5,436,134, and 5,582,977, and pending patent applications, and related foreign patents, for internal research and development use only in detecting, quantitating, and sizing macromolecules, in combination with microfluidics, where internal research and development use expressly excludes the use of this product for providing medical, diagnostic, or any other testing, analysis, or screening services, or providing clinical information or clinical analysis, in any event in return for compensation by an unrelated party.

Coomassie is a trademark of BASF Aktiengesellschaft.

Notice regarding Bio-Rad thermal cyclers and real-time systems.

Purchase of this instrument conveys a limited non-transferable immunity from suit for the purchaser's own internal research and development and for use in applied fields other than Human In Vitro Diagnostics under one or more of U.S. Patents Nos. 5,656,493, 5,333,675, 5,475,610 (claims 1, 44, 158, 160–163 and 167 only), and 6,703,236 (claims 1–7 only), or corresponding claims in their non-U.S. counterparts, owned by Applied Biosystems Corporation. No right is conveyed expressly, by implication or by estoppel under any other patent claim, such as claims to apparatus, reagents, kits, or methods such as 5' nuclease methods. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

Bio-Rad's real-time thermal cyclers are licensed real-time thermal cyclers under Applied Biosystems' United States Patent No. 6,814,934 B1 for use in research and for all other fields except the fields of human diagnostics and veterinary diagnostics.



**Bio-Rad
Laboratories, Inc.**

Life Science
Group

Web site www.bio-rad.com USA 800 4BIORAD Australia 02 9914 2800 Austria 01 877 89 01 Belgium 09 385 55 11 Brazil 55 21 3237 9400 Canada 905 712 2771 China 86 21 6426 0808 Czech Republic 420 241 430 532 Denmark 44 52 10 00 Finland 09 804 22 00 France 01 47 95 69 65 Germany 089 318 84 0 Greece 30 210 777 4396 Hong Kong 852 2789 3300 Hungary 36 1 455 8800 India 91 124 4029300/5013478 Israel 03 963 6050 Italy 39 02 216091 Japan 03 5811 6270 Korea 82 2 3473 4460 Mexico 55 5200 05 20 The Netherlands 0318 540666 New Zealand 64 9415 2280 Norway 23 38 41 30 Poland 48 22 331 99 99 Portugal 351 21 472 7700 Russia 7 095 721 14 04 Singapore 65 6415 3188 South Africa 27 0861 246 723 Spain 34 91 590 5200 Sweden 08 555 12700 Switzerland 061 717 95 55 Taiwan 886 2 2578 7189/2578 7241 United Kingdom 020 8328 2000